

# PIGS antibody



Catalog Number: 113897

## Product name

PIGS antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PIGS Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, HepG2 cells, mouse liver tissue, PC-3 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of PIGS recombinant protein (Accession Number: BC001319). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113897(PIGS Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:113897(PIGS Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)