

PIGT antibody

Catalog Number: 113898

Product name

PIGT antibody

Specificity

Human; other species not tested.

Antibody description

PIGT Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of PIGT recombinant protein (Accession Number: NM_015937). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113898(PIGT antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human brain using Catalog No:113898(PIGT antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:113898(PIGT antibody) at dilution of 1:100 (under 40x lens)