

# PIK3C2A antibody



Catalog Number: 113901

## Product name

PIK3C2A antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

## Specificity

Human, Mouse; other species not tested.


IHC: 1:20-1:200

IF: 1:50-1:500

## Antibody description

PIK3C2A Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human heart tissue, human placenta tissue. Positive IP detected in PC-3 cells. Positive WB detected in PC-3 cells, HeLa cells, human placenta tissue. Observed molecular weight by Western-blot: 191 kDa

## Validations



## Preparation

This antibody was obtained by immunization of PIK3C2A recombinant protein (Accession Number: XM\_017017924). Purification method: Antigen affinity purified.

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113901(PIK3C2A antibody) at dilution of 1:1500

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IP, IF

## Dilutions


Recommended Dilution:



Immunohistochemical of paraffin-embedded human heart using Catalog No:113901(PIK3C2A antibody) at dilution of 1:50 (under 10x lens)


# PIK3C2A antibody

Catalog Number: 113901



Immunohistochemical of paraffin-embedded human heart using Catalog No:113901(PIK3C2A antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-PIK3C2A (IP:Catalog No:113901, 4ug; Detection:Catalog No:113901 1:1000) with PC-3 cells lysate 3200ug.



Immunofluorescent analysis of HeLa cells using Catalog No:113901(PIK3C2A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

