

# PIK3IP1 antibody



Catalog Number: 113904

## Product name

PIK3IP1 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PIK3IP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse liver tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 37 kDa

## Preparation

This antibody was obtained by immunization of PIK3IP1 recombinant protein (Accession Number: NM\_052880). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113904(PIK3IP1 antibody) at dilution of 1:500



IP Result of anti-PIK3IP1 (IP:Catalog No:113904, 4ug; Detection:Catalog No:113904 1:500) with HeLa cells lysate 3200ug.