

PIK3IP1 antibody



Catalog Number: 113904

Product name

PIK3IP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PIK3IP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, mouse liver tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of PIK3IP1 recombinant protein (Accession Number: NM_052880). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113904(PIK3IP1 antibody) at dilution of 1:500



IP Result of anti-PIK3IP1 (IP:Catalog No:113904, 4ug; Detection:Catalog No:113904 1:500) with HeLa cells lysate 3200ug.