

PIP5K3 antibody



Catalog Number: 113916

Product name

Polyclonal

PIP5K3 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

PIP5K3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse brain tissue, mouse lung tissue. Observed molecular weight by Western-blot: 237kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of PIP5K3 recombinant protein (Accession Number: NM_152671). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113916(PIP5K3 antibody) at dilution of 1:300