

PKMYT1 antibody



Catalog Number: 113918

Product name

PKMYT1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PKMYT1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A375 cells, HeLa cells, HepG2 cells, HL-60 cells, mouse skin tissue. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004203). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113918(PKMYT1 antibody) at dilution of 1:500