

PLA2G1B antibody



Catalog Number: 113924

Product name

PLA2G1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLA2G1B Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 13-16 kDa

Preparation

This antibody was obtained by immunization of PLA2G1B recombinant protein (Accession Number: BC005386). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:113924(PLA2G1B Antibody) at dilution of 1:300