

PLAC1L antibody



Catalog Number: 113930

Product name

PLAC1L antibody

Specificity

Human; other species not tested.

Antibody description

PLAC1L Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Positive IHC detected in human ovary tissue, human placenta tissue. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of PLAC1L recombinant protein (Accession Number: NM_173801). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113930(PLAC1L Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human ovary tissue slide using Catalog No:113930(PLAC1L Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human ovary tissue slide using Catalog No:113930(PLAC1L Antibody) at dilution of 1:50

PLAC1L antibody



Catalog Number: 113930

(under 40x lens)