

PLAGL2 antibody



Catalog Number: 113937

Product name

PLAGL2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PLAGL2 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 50-55kd,70kd

Preparation

This antibody was obtained by immunization of PLAGL2 recombinant protein (Accession Number: XM_011528864). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:113937(PLAGL2 antibody) at dilution of 1:300