

# plasminogen antibody



Catalog Number: 113942

## Product name

plasminogen antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

plasminogen Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive IP detected in human plasma tissue. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 97 kDa

## Preparation

This antibody was obtained by immunization of plasminogen recombinant protein (Accession Number: NM\_000301). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:113942(PLG antibody) at dilution of 1:400



IP Result of anti-PLG (IP:Catalog No:113942, 2ug; Detection:Catalog No:113942 1:500) with human plasma tissue lysate 4000ug.




Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog

# plasminogen antibody

Catalog Number: 113942

---

No:113942(PLG Antibody) at dilution of 1:200  
(under 10x lens).



Immunohistochemistry of paraffin-embedded  
human liver tissue slide using Catalog  
No:113942(PLG Antibody) at dilution of 1:200  
(under 40x lens).