

PLCB3 antibody



Catalog Number: 113944

Product name

PLCB3 antibody

Specificity

Human; other species not tested.

Antibody description

PLCB3 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, A431 cells, COLO 320 cells, HEK-293T cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 150 kDa

Preparation

This antibody was obtained by immunization of PLCB3 recombinant protein (Accession Number: NM_000932). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113944(PLCB3 Antibody) at dilution of 1:600