

# PLCB3 antibody



Catalog Number: 113945

## Product name

Polyclonal

PLCB3 antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

PLCB3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293T cells, MCF-7 cells. Observed molecular weight by Western-blot: 139 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000932). Purification method: Antigen affinity purified.

## Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113945(PLCB3 Antibody) at dilution of 1:600