

# PLCB3 antibody

Catalog Number: 113945

---

**Product name**

Polyclonal

PLCB3 antibody

**Ig Type****Specificity**

Rabbit IgG

Human; other species not tested.

**Applications****Antibody description**

ELISA, WB

PLCB3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293T cells, MCF-7 cells. Observed molecular weight by Western-blot: 139 kDa


**Dilutions**

Recommended Dilution:

WB: 1:200-1:2000

**Preparation**

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000932). Purification method: Antigen affinity purified.

**Validations****Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113945(PLCB3 Antibody) at dilution of 1:600

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**