

# PLCD4 antibody

Catalog Number: 113948

---

**Product name**

PLCD4 antibody

**Specificity**

human, mouse, rat; other species not tested.

**Antibody description**

PLCD4 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer.

**Preparation**

This antibody was obtained by immunization of PLCD4 recombinant protein (Accession Number: NM\_032726). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, IHC

**Dilutions**

Recommended Dilution:

IHC: 1:20-1:200

**Validations**

Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113948(PLCD4 antibody) at dilution of 1:50 (under 10x lens)