

# PLD1 antibody

Catalog Number: 113951

## Product name

PLD1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PLD1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 120 kDa

## Preparation

This antibody was obtained by immunization of PLD1 recombinant protein (Accession Number: XM\_005247533). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113951(PLD1 antibody) at dilution of 1:500



IP Result of anti-PLD1 (IP:Catalog No:113951, 4ug; Detection:Catalog No:113951 1:500) with HeLa cells lysate 1200ug.



Immunofluorescent analysis of HeLa cells using Catalog No:113951(PLD1 Antibody) at dilution of

# PLD1 antibody



Catalog Number: 113951

---

1:50 and Alexa Fluor 488-conjugated AffiniPure

Goat Anti-Rabbit IgG(H+L)