

# PLEKHA1 antibody



Catalog Number: 113955

## Product name

PLEKHA1 antibody

## Immunogen

[Human PLEKHA1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PLEKHA1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in human placenta tissue, Jurkat cells. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of PLEKHA1 recombinant protein (Accession Number: XM\_017016481). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

## Validations

230kd →  
130kd →  
95kd →  
72kd →  
56kd →  
36kd →  
28kd →  
17kd →



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:113955(PLEKHA1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113955(PLEKHA1 antibody) at dilution of 1:100 (under 25x lens)