

# PLEKHA5 antibody



Catalog Number: 113957

## Product name

PLEKHA5 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PLEKHA5 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human heart tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 127kD,74kD

## Preparation

This antibody was obtained by immunization of PLEKHA5 recombinant protein (Accession Number: BC013133). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113957(PLEKHA5 antibody) at dilution of 1:500