

PLEKHO1 antibody



Catalog Number: 113964

Product name

PLEKHO1 antibody

Immunogen

[Human PLEKHO1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

PLEKHO1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse testis tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of PLEKHO1 recombinant protein (Accession Number: NM_016274). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113964(PLEKHO1 Antibody) at dilution of 1:300