

PLK3 antibody

Catalog Number: 113970

Product name

PLK3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLK3 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of PLK3 recombinant protein (Accession Number: NM_004073). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:10-1:100

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113970(PLK3 antibody) at dilution of 1:10 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113970(PLK3 Antibody) at dilution of 1:1000