

PNKD antibody

Catalog Number: 113973

Product name

PNKD antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PNKD Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, mouse brain tissue.

Preparation

This antibody was obtained by immunization of PNKD recombinant protein (Accession Number: NM_001077399). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:113973(PNKD antibody) at dilution of 1:100 (under 10x lens)