

PNMAL1 antibody



Catalog Number: 113978

Product name

PNMAL1 antibody

Specificity

human; other species not tested.

Antibody description

PNMAL1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A2780 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of PNMAL1 recombinant protein (Accession Number: BC026026). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113978(PNMAL1 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using PNMAL1 antibody Catalog No:113978 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).