

PNMT antibody

Catalog Number: 113979

Product name

PNMT antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PNMT Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, human heart tissue. Positive IP detected in PC-12 cells. Observed molecular weight by Western-blot: 31 kDa

Preparation

This antibody was obtained by immunization of PNMT recombinant protein (Accession Number: NM_002686). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:113979(PNMT antibody) at dilution of 1:400



IP Result of anti-PNMT (IP:Catalog No:113979, 4ug; Detection:Catalog No:113979 1:300) with PC-12 cells lysate 2400ug.