

PNRC1 antibody



Catalog Number: 113985

Product name

PNRC1 antibody

Specificity

human; other species not tested.

Antibody description

PNRC1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 35kd

Preparation

This antibody was obtained by immunization of PNRC1 recombinant protein (Accession Number: NM_006813). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:113985(PNRC1 antibody) at dilution of 1:500