

PODN antibody



Catalog Number: 113987

Product name

PODN antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PODN Rabbit Polyclonal antibody. Positive WB detected in A375 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of PODN recombinant protein (Accession Number: XM_011540677). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:113987(PODN Antibody) at dilution of 1:600