

# PLS3 antibody

Catalog Number: 114010

## Product name

PLS3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PLS3 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, COLO 320 cells, HeLa cells. Positive IHC detected in human oesophagus. Observed molecular weight by Western-blot: 63kd,38kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005032). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114010(PLS3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human oesophagus using Catalog No:114010(PLS3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human oesophagus using Catalog No:114010(PLS3 antibody) at dilution of 1:50 (under 40x lens)