

PLS3 antibody



Catalog Number: 114011

Product name

PLS3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PLS3 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 68 kDa

Preparation

This antibody was obtained by immunization of PLS3 recombinant protein (Accession Number: XM_047442170). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114011(PLS3 Antibody) at dilution of 1:600