

PMS1 antibody

Catalog Number: 114029

Product name

PMS1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PMS1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma.

Preparation

This antibody was obtained by immunization of PMS1 recombinant protein (Accession Number: BC008410). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114029(PMS1 antibody) at dilution of 1:100 (under 10x lens)