

PMS2 antibody



Catalog Number: 114030

Product name

PMS2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PMS2 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HeLa cells, HT-1080 cells, MCF7 cells, PC-3 cells, Raji cells. Observed molecular weight by Western-blot: 117 kDa, 100 kDa

Preparation

This antibody was obtained by immunization of PMS2 recombinant protein (Accession Number: NM_000535). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114030(PMS2 antibody) at dilution of 1:500