

POLR1D antibody



Catalog Number: 114037

Product name

Polyclonal

POLR1D antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

POLR1D Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, mouse lung tissue. Observed molecular weight by Western-blot: 20kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of POLR1D recombinant protein (Accession Number: XM_047430381). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114037(POLR1D antibody) at dilution of 1:500

Clonality