

# POLR1E antibody

Catalog Number: 114038

## Product name

POLR1E antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

POLR1E Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, human brain tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 47-53 kDa

## Preparation

This antibody was obtained by immunization of POLR1E recombinant protein (Accession Number: NM\_022490). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP, IHC

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114038(POLR1E antibody) at dilution of 1:500



IP Result of anti-POLR1E (IP:Catalog No:114038, 3ug; Detection:Catalog No:114038 1:500) with Jurkat cells lysate 1000ug.