

# POLR2B-Specific antibody



Catalog Number: 114039

## Product name

POLR2B-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

POLR2B-Specific Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in A431 cells, HeLa cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 140kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000938). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:200-1:1000

IHC: 1:20-1:200


## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114039(POLR2B-Specific antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:114039(POLR2B-Specific antibody) at dilution of 1:400 (under 10x



lens)

IP Result of anti-POLR2B-Specific (IP:Catalog No:114039, 4ug; Detection:Catalog No:114039 1:300) with HeLa cells lysate 1200ug.