

POLR2J2 antibody



Catalog Number: 114046

Product name

Polyclonal

POLR2J2 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

POLR2J2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Dilutions

Preparation

Recommended Dilution:

This antibody was obtained by immunization of POLR2J2 recombinant protein (Accession Number: NM_032958). Purification method: Antigen affinity purified.

WB: 1:500-1:5000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Validations



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114046(POLR2J2 Antibody) at dilution of 1:1000