

POLR3A antibody

Catalog Number: 114048

Product name

POLR3A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

POLR3A Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 140-150kd

Preparation

This antibody was obtained by immunization of POLR3A recombinant protein (Accession Number: NM_007055). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114048(POLR3A antibody) at dilution of 1:300