

# POLR3B antibody



Catalog Number: 114049

## Product name

POLR3B antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

POLR3B Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 128kd

## Preparation

This antibody was obtained by immunization of POLR3B recombinant protein (Accession Number: NM\_018082). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114049(POLR3B antibody) at dilution of 1:300