

POU3F4 antibody



Catalog Number: 114072

Product name

POU3F4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

POU3F4 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of POU3F4 recombinant protein (Accession Number: NM_000307). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114072(POU3F4 Antibody) at dilution of 1:300