

# POU3F4 antibody



Catalog Number: 114072

## Product name

POU3F4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

POU3F4 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 40-45 kDa

## Preparation

This antibody was obtained by immunization of POU3F4 recombinant protein (Accession Number: NM\_000307). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114072(POU3F4 Antibody) at dilution of 1:300