

POU4F2 antibody



Catalog Number: 114073

Product name

POU4F2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

POU4F2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, mouse eye tissue. Observed molecular weight by Western-blot: 40-44kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114073(POU4F2 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells using Catalog No:114073(POU4F2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)