

POU4F3 antibody



Catalog Number: 114074

Product name

POU4F3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

POU4F3 Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells, HL-60 cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of POU4F3 recombinant protein (Accession Number: NM_002700). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:114074(POU4F3 Antibody) at dilution of 1:1000