

# PPP1R9A antibody



Catalog Number: 114093

## Product name

PPP1R9A antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PPP1R9A Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 180kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001166160). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114093(PPP1R9A antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using PPP1R9A antibody Catalog No:114093 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).