

PPP2R1A antibody



Catalog Number: 114097

Product name

PPP2R1A antibody

Immunogen

[Human PPP2R1A Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPP2R1A Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A431 cells, HeLa cells, mouse kidney tissue. Positive IHC detected in human kidney tissue, human ovary tumor tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of PPP2R1A recombinant protein (Accession Number: NM_014225). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114097(PPP2R1A antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114097(PPP2R1A antibody) at dilution of 1:50 (under 10x lens)



PPP2R1A antibody



Catalog Number: 114097

Immunohistochemical of paraffin-embedded

human kidney using Catalog No:114097(PPP2R1A antibody) at dilution of 1:50 (under 40x lens)