

# PPP2R1B antibody



Catalog Number: 114098

---

## Product name

PPP2R1B antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PPP2R1B Rabbit Polyclonal antibody. Positive WB detected in A2780 cells, human kidney tissue, mouse brain tissue. Observed molecular weight by Western-blot: 74kd

## Preparation

This antibody was obtained by immunization of PPP2R1B recombinant protein (Accession Number: NM\_181699). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:114098(PPP2R1B antibody) at dilution of 1:500