

# PPP2R2A antibody

Catalog Number: 114099

## Product name

PPP2R2A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPP2R2A Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 50 kDa, 62 kDa

## Preparation

This antibody was obtained by immunization of PPP2R2A recombinant protein (Accession Number: NM\_002717). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114099(PPP2R2A Antibody) at dilution of 1:300



IP Result of anti-PPP2R2A (IP:Catalog No:114099, 4ug; Detection:Catalog No:114099 1:500) with mouse brain tissue lysate 4000ug.