

PPP3CA-specific antibody



Catalog Number: 114109

Product name

PPP3CA-specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPP3CA-specific Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse pancreas tissue. Positive IHC detected in human heart tissue, human brain tissue. Observed molecular weight by Western-blot: 59 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000944). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114109(PPP3CA antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human heart using Catalog No:114109(PPP3CA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:114109(PPP3CA antibody) at dilution of 1:50 (under 40x lens)