

# PPP3CC antibody



Catalog Number: 114112

## Product name

PPP3CC antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPP3CC Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue, MCF7 cells, mouse brain tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of PPP3CC recombinant protein (Accession Number: NM\_005605). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114112(PPP3CC antibody) at dilution of 1:500