

PPP3CC-specific antibody



Catalog Number: 114113

Product name

PPP3CC-specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PPP3CC-specific Rabbit Polyclonal antibody.
Positive WB detected in fetal human brain tissue, human brain tissue, human testis tissue. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005605).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114113(PPP3CC Antibody) at dilution of 1:1000