

PPEF1 antibody



Catalog Number: 114119

Product name

PPEF1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPEF1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human brain tissue. Observed molecular weight by Western-blot: 76 kDa

Preparation

This antibody was obtained by immunization of PPEF1 recombinant protein (Accession Number: XM_005274553). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114119(PPEF1 antibody) at dilution of 1:300