

# PPFIBP2 antibody



Catalog Number: 114122

---

## Product name

PPFIBP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PPFIBP2 Rabbit Polyclonal antibody. Positive IF detected in A431 cells.

## Preparation

This antibody was obtained by immunization of PPFIBP2 recombinant protein (Accession Number: NM\_001351862). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of A431 cells using Catalog No:114122(PPFIBP2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)