

PPIL4 antibody



Catalog Number: 114130

Product name

PPIL4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPIL4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 58-65kd

Preparation

This antibody was obtained by immunization of PPIL4 recombinant protein (Accession Number: BC016984). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114130(PPIL4 antibody) at dilution of 1:1000