

PPM1F antibody



Catalog Number: 114134

Product name

PPM1F antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPM1F Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, HepG2 cells, human liver tissue, K-562 cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of PPM1F recombinant protein (Accession Number: NM_014634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114134(PPM1F antibody) at dilution of 1:500